

Product datasheet for **TP726445**

BPGM Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Bisphosphoglycerate mutase/BPGM
Species:	Human
Expression cDNA Clone or AA Sequence:	G108-K259
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Bisphosphoglycerate mutase/BPGM is produced by E. coli. The target gene encoding G108-K259 is expressed with a FLAG tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	669
UniProt ID:	P07738



[View online »](#)