

Product datasheet for **TP726447**

CENPE Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Centromere-associated protein E/CENPE
Species:	Human
Expression cDNA Clone or AA Sequence:	K2400-R2500
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Centromere-associated protein E/CENPE is produced by E. coli. The target gene encoding K2400-R2500 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1062
UniProt ID:	Q02224



[View online »](#)