

## Product datasheet for **TP726448**

### **C11orf54 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Ester hydrolase C11orf54/C11orf54
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	A2-G100
<b>Tag:</b>	C-FLAG
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human Ester hydrolase C11orf54/C11orf54is produced by E. coli. The target gene encodingA2-G100 is expressed with a FLAG tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	28970
<b>UniProt ID:</b>	<u><a href="#">Q9H0W9</a></u>



[View online »](#)