

Product datasheet for **TP726449**

Proteasome 20S beta 6 (PSMB6) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Proteasome subunit beta type-6/PSMB6
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-A239
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Proteasome subunit beta type-6/PSMB6is produced by E. coli. The target gene encodingA2-A239 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5694
UniProt ID:	P28072



[View online »](#)