

## Product datasheet for **TP726450**

### ARL2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ADP-ribosylation factor-like protein 2/ARL2
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-D184
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ADP-ribosylation factor-like protein 2/ARL2is produced by E. coli. The target gene encodingG2-D184 is expressed with a FLAG tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	402
UniProt ID:	<a href="#">P36404</a>
Synonyms:	ADP-ribosylation factor-like protein 2; ARL2



[View online »](#)