

Product datasheet for **TP726450**

ARL2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ADP-ribosylation factor-like protein 2/ARL2
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-D184
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ADP-ribosylation factor-like protein 2/ARL2is produced by E. coli. The target gene encodingG2-D184 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	402
UniProt ID:	P36404



[View online »](#)