

Product datasheet for **TP726451**

Acidic Calponin (CNN3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Calponin-3/CNN3
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-Y329
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Calponin-3/CNN3is produced by E. coli. The target gene encodingT2-Y329 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1266
UniProt ID:	<u>Q15417</u>



[View online »](#)