

## Product datasheet for **TP726455**

### NUDCD2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NudC domain-containing protein 2/NUDCD2
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-K157
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human NudC domain-containing protein 2/NUDCD2is produced by E. coli. The target gene encodingS2-K157 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	134492
UniProt ID:	<a href="#">Q8WVJ2</a>



[View online »](#)