

Product datasheet for **TP726456**

Malectin (MLEC) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Malectin/MLEC
Species:	Human
Expression cDNA Clone or AA Sequence:	E106-S269
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Malectin/MLEC is produced by E. coli. The target gene encoding E106-S269 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	9761
UniProt ID:	Q14165



[View online »](#)