

Product datasheet for **TP726458**

VPS29 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human VPS29 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	N68-P182
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human VPS29 Protein is produced by E. coli. The target gene encoding N68-P182 is expressed with a FLAG tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51699
UniProt ID:	Q9UBQ0



[View online »](#)