

Product datasheet for **TP726459**

CDKL2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cyclin-dependent kinase-like 2/CDKL2
Species:	Human
Expression cDNA Clone or AA Sequence:	N220-K347
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Cyclin-dependent kinase-like 2/CDKL2 is produced by E. coli. The target gene encodingN220-K347 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8999
UniProt ID:	<u>Q92772</u>



[View online »](#)