

Product datasheet for **TP726465**

LRRFIP2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LRRFIP2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	H524-Y661
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human LRRFIP2 Protein is produced by E. coli. The target gene encoding H524-Y661 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	9209
UniProt ID:	Q9Y608



[View online »](#)