

Product datasheet for **TP726466**

OAZ1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ornithine decarboxylase antizyme 1/OAZ1
Species:	Human
Expression cDNA Clone or AA Sequence:	F95-E219
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ornithine decarboxylase antizyme 1/OAZ1 is produced by E. coli. The target gene encoding F95-E219 is expressed with a FLAG tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4946
UniProt ID:	P54368
Synonyms:	Ornithine decarboxylase antizyme 1; OAZ1; AZ1; ODC-Az; OAZ



[View online »](#)