

Product datasheet for **TP726467**

MTM1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Myotubularin/MTM1
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-T100
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Myotubularin/MTM1 is produced by E. coli. The target gene encoding A2-T100 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4534
UniProt ID:	Q13496



[View online »](#)