

Product datasheet for **TP726469**

Factor VIII (F8) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Coagulation factor VIII/F8 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | I1519-F1718 |
| Tag: | C-FLAG |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | Recombinant Human Coagulation factor VIII/F8is produced by E. coli. The target gene encodingI1519-F1718 is expressed with a FLAG tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 2157 |
| UniProt ID: | <u>P00451</u> |



[View online »](#)