

Product datasheet for **TP726471**

BMP6 Human Recombinant Protein

Product data:

| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Bone morphogenetic protein 6/BMP6 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | V370-H513 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | Recombinant Human Bone morphogenetic protein 6/BMP6is produced by E. coli. The target gene encodingV370-H513 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 654 |
| UniProt ID: | <u>P22004</u> |



[View online »](#)