

## Product datasheet for **TP726473**

### Caspase 3 (CASP3) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Caspase-3/CASP3
Species:	Human
Expression cDNA Clone or AA Sequence:	S29-K210
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Caspase-3/CASP3is produced by E. coli. The target gene encodingS29-K210 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	836
UniProt ID:	<a href="#">P42574</a>



[View online »](#)