

Product datasheet for **TP726474**

PAFAH (PLA2G7) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Platelet-activating factor acetylhydrolase/PLA2G7
Species:	Human
Expression cDNA Clone or AA Sequence:	K210-N441
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Platelet-activating factor acetylhydrolase/PLA2G7is produced by E. coli. The target gene encodingK210-N441 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7941
UniProt ID:	<u>Q13093</u>



[View online »](#)