

Product datasheet for **TP726475**

IL17D Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Interleukin-17D/IL17D
Species:	Human
Expression cDNA Clone or AA Sequence:	A18-P202
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Interleukin-17D/IL17Dis produced by E. coli. The target gene encodingA18-P202 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	53342
UniProt ID:	Q8TAD2
Synonyms:	Interleukin-17D; IL17D; Interleukin-27 (IL-27); UNQ3096/PRO21175
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)