

Product datasheet for **TP726488**

IL32 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Interleukin-32/IL32
Species:	Human
Expression cDNA Clone or AA Sequence:	C2-K131
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Interleukin-32/IL32is produced by E. coli. The target gene encodingC2-K131 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9235
UniProt ID:	P24001
Synonyms:	Interleukin-32; Natural killer cells protein 4; Tumor necrosis factor alpha-inducing factor; IL32; NK4; TAIF
Protein Families:	Secreted Protein



[View online »](#)