

## Product datasheet for **TP726490**

### Prolactin (PRL) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Prolactin/PRL
Species:	Human
Expression cDNA Clone or AA Sequence:	L60-C227
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Prolactin/PRL is produced by E. coli. The target gene encoding L60-C227 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5617
UniProt ID:	<a href="#">P01236</a>



[View online »](#)