

## Product datasheet for **TP726492**

### XIAP Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human E3 ubiquitin-protein ligase XIAP
Species:	Human
Expression cDNA Clone or AA Sequence:	E350-S497
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human E3 ubiquitin-protein ligase XIAP is produced by E. coli. The target gene encoding E350-S497 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	331
UniProt ID:	<a href="#">P98170</a>



[View online »](#)