

Product datasheet for **TP726493**

S adenosylhomocysteine hydrolase (AHCY) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Adenosylhomocysteinase/AHCY
Species:	Human
Expression cDNA Clone or AA Sequence:	V180-A340
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Adenosylhomocysteinase/AHCY is produced by E. coli. The target gene encoding V180-A340 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	191
UniProt ID:	P23526



[View online »](#)