

Product datasheet for **TP726494**

Calcineurin B (PPP3R1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Calcineurin subunit B type 1/PPP3R1
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-V170
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Calcineurin subunit B type 1/PPP3R1 is produced by E. coli. The target gene encoding G2-V170 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5534
UniProt ID:	P63098



[View online »](#)