

## Product datasheet for **TP726498**

### FRMPD2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FERM and PDZ domain-containing protein 2/FRMPD2
Species:	Human
Expression cDNA Clone or AA Sequence:	R193-P311
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human FERM and PDZ domain-containing protein 2/FRMPD2is produced by E. coli. The target gene encodingR193-P311 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	143162
UniProt ID:	<a href="#">Q68DX3</a>



[View online »](#)