

## Product datasheet for **TP726500**

### Endonuclease V (ENDO V) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Endonuclease V/ENDO V
Species:	Human
Expression cDNA Clone or AA Sequence:	V153-C282
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Endonuclease V/ENDO V is produced by E. coli. The target gene encoding V153-C282 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	284131
UniProt ID:	<a href="#">Q8N8Q3</a>
Synonyms:	Endonuclease V; hEndoV; Inosine-specific endoribonuclease; ENDO V
Protein Families:	Druggable Genome



[View online »](#)