

Product datasheet for **TP726503**

KLHDC10 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Kelch domain-containing protein 10/KLHDC10
Species:	Human
Expression cDNA Clone or AA Sequence:	G6-A155
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Kelch domain-containing protein 10/KLHDC10 is produced by E. coli. The target gene encodingG6-A155 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23008
UniProt ID:	Q6PID8
Synonyms:	Kelch domain-containing protein 10; KLHDC10; KIAA0265



[View online »](#)