

Product datasheet for **TP726504**

ASB8 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ankyrin repeat and SOCS box protein 8/ASB8
Species:	Human
Expression cDNA Clone or AA Sequence:	Q13-P121
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ankyrin repeat and SOCS box protein 8/ASB8is produced by E. coli. The target gene encodingQ13-P121 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	140461
UniProt ID:	Q9H765



[View online »](#)