

## Product datasheet for **TP726504**

### ASB8 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ankyrin repeat and SOCS box protein 8/ASB8
Species:	Human
Expression cDNA Clone or AA Sequence:	Q13-P121
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ankyrin repeat and SOCS box protein 8/ASB8 is produced by E. coli. The target gene encoding Q13-P121 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	140461
UniProt ID:	<a href="#">Q9H765</a>
Synonyms:	Ankyrin repeat and SOCS box protein 8; ASB8; PP14212
Protein Families:	Druggable Genome



[View online »](#)