

Product datasheet for **TP726505**

ACTL9 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin-like protein 9/ACTL9
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-S172
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Actin-like protein 9/ACTL9 is produced by E. coli. The target gene encodingD2-S172 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	284382
UniProt ID:	Q8TC94



[View online »](#)