

Product datasheet for **TP726510**

CCDC42B (CFAP73) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cilia- and flagella-associated protein 73/CFAP73
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-R100
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Cilia- and flagella-associated protein 73/CFAP73 is produced by E. coli. The target gene encoding A2-R100 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	387885
UniProt ID:	A6NFT4
Synonyms:	Cilia- and flagella-associated protein 73; CFAP73; Coiled-coil domain-containing protein 42B; CCDC42B



[View online »](#)