

Product datasheet for **TP726511**

TBC1D8B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TBC1 domain family member 8B/TBC1D8B
Species:	Human
Expression cDNA Clone or AA Sequence:	L936-A1084
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human TBC1 domain family member 8B/TBC1D8Bis produced by E. coli. The target gene encodingL936-A1084 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54885
UniProt ID:	Q0IIM8
Synonyms:	TBC1 domain family member 8B; TBC1D8B



[View online »](#)