

## Product datasheet for **TP726516**

### CDKL4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cyclin-dependent kinase-like 4/CDKL4
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-K104
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Cyclin-dependent kinase-like 4/CDKL4is produced by E. coli. The target gene encodingE2-K104 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	344387
UniProt ID:	<a href="#">Q5MAI5</a>
Synonyms:	Cyclin-dependent kinase-like 4; CDKL4
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)