

## Product datasheet for **TP726517**

### **LDHAL6A Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human L-lactate dehydrogenase A-like 6A/LDHAL6A
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	P210-L332
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human L-lactate dehydrogenase A-like 6A/LDHAL6A is produced by E. coli. The target gene encoding P210-L332 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	160287
<b>UniProt ID:</b>	<a href="#">Q6ZMR3</a>



[View online »](#)