

## Product datasheet for **TP726520**

### EXOC3L1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Exocyst complex component 3-like protein/EXOC3L1
Species:	Human
Expression cDNA Clone or AA Sequence:	G352-L488
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Exocyst complex component 3-like protein/EXOC3L1 is produced by E. coli. The target gene encoding G352-L488 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	283849
UniProt ID:	<a href="#">Q86V11</a>
Synonyms:	Exocyst complex component 3-like protein; EXOC3L1; Protein Jiangli; EXOC3L



[View online »](#)