

Product datasheet for **TP726522**

MTURN Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Maturin/MTURN
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-Q131
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Maturin/MTURN is produced by E. coli. The target gene encoding D2-Q131 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	222166
UniProt ID:	Q8N3F0
Synonyms:	Maturin; Maturin neural progenitor differentiation regulator protein homolog; Protein Ells1; MTURN; C7orf41



[View online »](#)