

Product datasheet for **TP726523**

DOC2A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human DOC2A Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A137-Q244
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human DOC2A Protein is produced by E. coli. The target gene encoding A137-Q244 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8448
UniProt ID:	Q14183
Synonyms:	Double C2-like domain-containing protein alpha; DOC2A; Doc2; Doc2-alpha
Protein Families:	Secreted Protein



[View online »](#)