

Product datasheet for **TP726527**

Speedy protein C (SPDYC) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Speedy protein C/SPDYC
Species:	Human
Expression cDNA Clone or AA Sequence:	D120-H293
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Speedy protein C/SPDYC is produced by E. coli. The target gene encodingD120-H293 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	387778
UniProt ID:	Q5MJ68
Synonyms:	Speedy protein C; Rapid inducer of G2/M progression in oocytes C (RINGO C; hSpy/Ringo C); SPDYC
Protein Pathways:	Oocyte meiosis, Progesterone-mediated oocyte maturation



[View online »](#)