

## Product datasheet for **TP726528**

### **APRIL (TNFSF13) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human TNFSF13 Protein
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	A105-L250
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human TNFSF13 Protein is produced by E. coli. The target gene encoding A105-L250 is expressed with a His tag at the C-terminus.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	8741
<b>UniProt ID:</b>	<a href="#">O75888</a>



[View online »](#)