

## Product datasheet for **TP726529**

### TRAF2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TNF receptor-associated factor 2/TRAF2
Species:	Human
Expression cDNA Clone or AA Sequence:	V40-A220
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human TNF receptor-associated factor 2/TRAF2is produced by E. coli. The target gene encodingV40-A220 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7186
UniProt ID:	<a href="#">Q12933</a>



[View online »](#)