

Product datasheet for **TP726531**

RASL10A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras-like protein family member 10A/RASL10A
Species:	Human
Expression cDNA Clone or AA Sequence:	S100-M203
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ras-like protein family member 10A/RASL10A is produced by E. coli. The target gene encoding S100-M203 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10633
UniProt ID:	Q92737
Synonyms:	Ras-like protein family member 10A; RASL10A; Ras-like protein RRP22; Ras-related protein on chromosome 22; RRP22



[View online »](#)