

Product datasheet for **TP726535**

EMX2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Homeobox protein EMX2
Species:	Human
Expression cDNA Clone or AA Sequence:	Q112-D252
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Homeobox protein EMX2is produced by E. coli. The target gene encodingQ112-D252 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	2018
UniProt ID:	Q04743



[View online »](#)