

Product datasheet for **TP726538**

NDUF3 (NDUFAF3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NDUFAF3 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A37-C150
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human NDUFAF3 Protein is produced by E. coli. The target gene encoding A37-C150 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	25915
UniProt ID:	Q9BU61



[View online »](#)