

Product datasheet for **TP726539**

TAF9B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TAF9B Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-M251
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human TAF9B Protein is produced by E. coli. The target gene encoding E2-M251 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	51616
UniProt ID:	<u>Q9HBM6</u>



[View online »](#)