

Product datasheet for **TP726546**

FEM1A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein fem-1 homolog A/FEM1A
Species:	Human
Expression cDNA Clone or AA Sequence:	R100-A300
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Protein fem-1 homolog A/FEM1A is produced by E. coli. The target gene encoding R100-A300 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55527
UniProt ID:	Q9BSK4
Synonyms:	Protein fem-1 homolog A; FEM1a; FEM1-alpha; Prostaglandin E receptor 4-associated protein; FEM1A; EPRAP
Protein Families:	Druggable Genome



[View online »](#)