

## Product datasheet for **TP726548**

### Neurogenin3 (NEUROG3) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Neurogenin-3/NEUROG3
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-L214
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Neurogenin-3/NEUROG3is produced by E. coli. The target gene encodingT2-L214 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	50674
UniProt ID:	<a href="#">Q9Y4Z2</a>



[View online »](#)