

Product datasheet for **TP726549**

SMUG1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SMUG1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	T116-K270
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human SMUG1 Protein is produced by E. coli. The target gene encoding T116-K270 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23583
UniProt ID:	Q53HV7
Synonyms:	Single-strand selective monofunctional uracil DNA glycosylase; SMUG1
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair



[View online »](#)