

Product datasheet for **TP726557**

TLX2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human T-cell leukemia homeobox protein 2/TLX2
Species:	Human
Expression cDNA Clone or AA Sequence:	W120-V284
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human T-cell leukemia homeobox protein 2/TLX2 is produced by E. coli. The target gene encoding W120-V284 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3196
UniProt ID:	O43763
Synonyms:	T-cell leukemia homeobox protein 2; Homeobox protein Hox-11L1; Neural crest homeobox protein; TLX2; HOX11L1; NCX
Protein Families:	Transcription Factors



[View online »](#)