

## Product datasheet for **TP726558**

### USP26 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 26/USP26
Species:	Human
Expression cDNA Clone or AA Sequence:	T664-E913
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 26/USP26is produced by E. coli. The target gene encodingT664-E913 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	83844
UniProt ID:	<a href="#">Q9BXU7</a>



[View online »](#)