

Product datasheet for **TP726560**

HMX2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Homeobox protein HMX2
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-A215
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Homeobox protein HMX2is produced by E. coli. The target gene encodingG2-A215 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3167
UniProt ID:	A2RU54
Synonyms:	Homeobox protein HMX2; HMX2; Homeobox protein H6 family member 2



[View online »](#)